

1

Non-viral Delivery of Prime Editing (PE): Analytical Tools for Testing and Characterizing PE RNAs-LNP and Case Studies

> Minjia Wang, PhD | Scientist, Prime Medicine On behalf of the team at Prime Medicine





• Minjia Wang declares he is a current employee of Prime Medicine, Inc. and owns equity in Prime Medicine







- 1. Gene therapy using LNP delivery of Prime Editing (PE)
- 2. LNP characteristics and attributes
- 3. Analytical methods available for PE LNP drug products
- 4. Analytical assays supporting LNP drug products developed by Prime:
  - a. LNP lipid component ratio analysis
  - b. LNP payload integrity and component ratio analysis



#### Lipid Nanoparticles (LNP) – A Versatile Delivery Method



- LNP with ionizable lipid could effectively encapsulate nucleic acids with increased stability, enable CRISPR-based gene editing both *ex vivo* and *in vivo*
- Major advantages: tunable, tissue/cell targeting potential, versatile, scalable, and no genome integration compared with viral gene delivery



### Prime Editing: A Precise Gene Editing Technology That Is Programmable for Both Search and Replace



5

• Prime Editing can deliver a one-time curative therapy that has the potential to address a large percentage of genetic variants associated with disease



pegRNA: prime editing guide RNA; RT: reverse transcriptase; Cas: CRISPR associated protein;

Anzalone, et al (David R. Liu). Search-and-replace genome editing without double-strand breaks or donor DNA. Nature, 2019.

Therapeutic Approach: LNP-mediated Targeting Delivery of Prime Editing to Liver



- Prime Medicine is developing a universal liver targeting RNA-LNP delivery platform that has the potential to expedite development of Prime Editing therapeutics
- Prime Editing can be used to new liver indication by just swapping the guide RNAs in the modular universal LNP and leveraging certain CMC, off-target and PharmTox packages across programs
- LNP co-delivers 3 Prime Editing components: pegRNA, ngRNA, and PE encoding mRNA



# Current Prime Medicine Pipeline Aims to Leverages the Versatility and Breadth of Prime Editing



• Universal liver targeting RNA-LNP delivery platform to be used for different liver indications

STRATEGIC CATEGORY	TARGET TISSUE	INDICATION	DELIVERY	DISCOVERY	LEAD OPTIMIZATION	IND-ENABLING	Phase 1/2
STRATEGIC CATEGORY IMMEDIATE	BLOOD	Chronic Granulomatous Disease	ex vivo				
		Fanconi Anemia	ex vivo				
	LIVER	Wilson's Disease	LNP				
		Glycogen Storage Disease 1b	LNP				
	EYE	Retinitis Pigmentosa/Rhodopsin	AAV				
		Retinitis Pigmentosa/Usher Syndrome	AAV				
	EAR	Usher Syndrome Type 3	AAV				
		Non-Syndromic Hearing Loss – GJB2	AAV				
	NEURO- MUSCULAR	Friedreich's Ataxia	AAV				
		Myotonic Dystrophy Type 1	viral/non-viral				
DIFFERENTIATION: REPEAT EXPANSION DISEASES		Amyotrophic Lateral Sclerosis	viral/non-viral				
		Oculopharyngeal Muscular Dystrophy	LNP				
		Fragile X Syndrome	viral/non-viral				
		Huntington's Disease	TBD				
	EYE	Fuchs' Endothelial Corneal Dystrophy	viral/non-viral				
DIFFERENTIATION: OTHER	MUSCLE	Duchenne Muscular Dystrophy	AAV				
	LUNG	Cystic Fibrosis	LNP				

BLOOD



# Analytical Assays Used for LNP Delivered PE Drug Product to the Liver medicine



- Battery of assays required for release, stability and characterization
- Many assays for gRNA, mRNA, and DP are developed as platform assays to support multiple programs
- Product specific potency assays are being developed early, where possible, to streamline CMC and clinical data correlation

#### Lipid Nanoparticles (LNP) – Characteristics and Analytical Tools

• Important attributes allow for the characterization of LNP based therapeutics

#### LNP specific attributes:

- Morphology (CryoTEM, SAXS, SANS, etc.)
- Biophysical: size and PDI (DLS, FFF/SEC-MALS, nano flow cytometer, etc.), surface charge, pKa
- Payload encapsulation efficiency
- Lipid (LC-MS/ELSD/CAD) and cargo (IP-RP-LC-UV, CGE/FA, NGS/ddPCR) ID, integrity, quantity

Other attributes:

- Safety: endotoxin, sterility, etc.
- Process related impurities
- Potency/functional activity
- Compendial: appearance, osmolarity, pH, etc.
- Stability

PDI: Polydispersity; FA: Fragment analyzer; CAD: Charged aerosol detector ELSD: Evaporative light scattering detector; NGS: Next gen sequencing

TEM: Transmission electron microscopy; DLS: Dynamic light scattering; SAXS: Small-angle X-ray scattering SANS: Small-angle neutron scattering; SEC-MALS: Size exclusion chromatography-multi-angle light scattering;

9

LNP gene delivery LNP analysis methods

Lipid analysis

sis RNA integrity analysis

RNA ratio analysis





prime\_ medicine

#### LNP Composition/Lipid Ratio Analysis Using RP-UPLC-ELSD/MS



- Important parameters: diversity in lipid components, chromatography mode, sample preparation, etc.
- Fast (5 min run), sensitive, and robust reversed-phased (RP) LC method using both ELSD and MS (Tof).
- Well resolved peaks without matrix interference.



#### LNP Composition/Lipid Ratio Analysis Using RP-UPLC-ELSD/MS



• Optimized LC conditions that can easily accommodate different lipid components, ratios, and concentrations in different LNP samples



#### LNP Composition/Lipid Ratio Analysis Using RP-UPLC-ELSD/MS



• Accurate quantification using both ELSD and MS detector for LNP lipid ratio analysis with high consistency to expected values



Representative lipid ratio analysis in disrupted LNP sample

Component	Ionizable lipid	Helper lipid	Sterol	PEG- lipid
Normalized expected lipid concentration (%)	100	100	100	100
Actual molar concentration as % of expected value	99.2	104.6	98.3	109.0

**RNA** ratio analysis

RNA integrity analysis

#### Impurity/Purity Analysis in Lipid Raw Material Using LC-MS

- prime\_ medicine
- Successful separation, structural identification, formation mechanism, and relative quantification of impurities in raw material (ionizable lipid) before and after purification



13

Lipid analysis

RNAs Payload Analysis: Impurity Analysis with CGE and RNA Quantification/Ratio Analysis with IP-RP-HPLC-UV



• Optimized workflow for sample preparation and payload analysis



Capillary gel electrophoresis (CGE)

Ion paring-reversed phase-HPLC-UV (IP-RP-HPLC-UV) High Resolution and Sensitivity Integrity Analysis of the Prime Editing Components in LNP with Capillary Gel Electrophoresis



- Well resolved gRNA components (pegRNA and ngRNA) with  $\sim 20$  n.t. difference by CGE ٠
- Streamlined LNP disruption method, fluorescence dye selection, and optimized CGE run conditions ٠ enable high sensitivity integrity analysis in different LNP samples



15

#### RNA Components Ratio/Quantification Analysis Using IP-RP-HPLC

prime\_ medicine

- Optimized IP-RP-LC conditions can resolve different gRNA components with ~ 20 n.t. difference
- Calibration curves with good linearity for accurate RNA quantification in different LNP samples



as % of expected value

Lipid analysis

**RNA** ratio analysis

#### Summary

- Prime Editing is considered a third generation gene editing technology that has the potential to address a large percentage of genetic variants associated with disease
- Prime Medicine is developing a universal LNP-RNA platform to deliver Prime Editors to treat liver and metabolic diseases
- Prime Medicine has developed multiple platform analytical assays that can support different LNP programs
- Three platform analytical assays developed internally by Prime Medicine were shared to show the lipid and payload analysis of LNP drug products

## **THANK YOU!**



